

S. Lefèvre<sup>1</sup>, Y. Ruano<sup>2</sup>, J.L. Rodríguez Peralto<sup>2</sup>, F. Rojo Todo<sup>3</sup>, N. Carvajal García<sup>3</sup>, I. Abdulkader Nallib<sup>4</sup>, M. Sánchez Ares<sup>4</sup>, J. Hernández Losa<sup>5</sup>, M. Sese Faustino<sup>5</sup>, K. Ashour Garrido<sup>1</sup>, M. Liénard<sup>1</sup>, J. Hong<sup>1</sup>, M. Elias<sup>1</sup>, A.B. Rodrigo Martínez<sup>1</sup>, A. Felipe Benavent<sup>1</sup>  
<sup>1</sup> OncoDNA SA, Rue Louis Breguet 1, Gosselies, Belgium. <sup>2</sup> University Hospital October 12, Servicio de Anatomía Patológica, Laboratorio de Patología Molecular, Madrid, Spain. <sup>3</sup> University Hospital Fundación Jiménez Díaz, Servicio de Patología Molecular, Madrid, Spain. <sup>4</sup> Santiago Clinic Hospital CHUS, Servicio de Anatomía Patológica, Santiago de Compostela, Spain. <sup>5</sup> University Hospital Vall d'Hebron, Servicio de Anatomía Patológica - Laboratorio de Patología Molecular, Barcelona, Spain

## Introduction :

- Homologous Recombination Repair (HRR) pathway restores DNA double-strand breaks and deficiency in this pathway is called **Homologous recombination deficiency (HRD)** which leads to genomic instability
- HRD biomarker can **predict response** to specific therapies, such as PARP inhibitors, and helps clinician to determine **the best treatment option** for the patients
- OncoDEEP® kit is a Comprehensive Genomic Panel allowing the analysis of both **DNA and RNA** combined with **Bio-IT analysis and clinical reporting**.
- Complete tumor characterization **through 638 genes, genomic signatures and fusions (Figure 1)** that helps **reducing the costs of testing and deliver faster results** in the selection of appropriate cancer treatment options.




Figure 1: Alterations and biomarkers analyzed by the OncoDEEP panel at the DNA level.

## Material and Methods:

- 76 ovarian carcinoma FFPE samples were selected from 4 Spanish Hospitals according to study criteria (Figure 2) and OncoDEEP kit was performed in the laboratory facilities of those 4 centers.




Figure 2: Selection criteria of 76 samples for the HRD OncoDEEP analytical validation and distribution of those samples in the 4 spanish centers

- The OncoDEEP kit allowed the full characterization, from the DNA extraction to the final report, of the samples in less than 5 working days (Figure 3). After extraction, libraries were constructed (3h), enriched (hands on time: 4h) based on Twist Biosciences Technology and sequenced (20h) on an Illumina NextSeq 500 or NextSeq 2000, depending on the center. Finally, FastQ files were uploaded and analyzed through OncoDNA dedicated BioIT pipeline.




Figure 3: Full OncoDEEP workflow from the wetlab part (DNA extraction to sequencing) to FASTQ files upload and secondary/tertiary analysis.

- HRD results were obtained combining the BRCA1/2 status and the Genomic Scar (GS) (Figure 4). The three biomarkers Loss Of Heterozygosity (LOH), Allelic Disparity on Telomere (ADT) and Large-scale Rearrangements (LR) were used for GS calculation. Then, those results were compared with the one obtained with Myriad MyChoice DX. An overall concordance >80% and BRCA1/2 status concordance >90% were to be reached.




Figure 4: HRD determination by OncoDEEP combining both BRCA1/2 status and GS and results criteria.

**Conclusion:** OncoDEEP® kit is an end-to-end solution allowing the detection of variants in 638 genes and the calculation of genomic signatures such as Tumor Mutational Burden (TMB), Microsatellite Instability (MSI) and HRD. OncoDEEP® HRD analysis showed an overall concordance of **88.8%** with a sensitivity and specificity of **90%** and **88%**, respectively, compared to the gold standard and is therefore applicable to the identification of ovarian cancer patients eligible for PARPi therapy.

## Results

- From the 76 samples (Table 1), 5 were rejected due to low coverage and/or uniformity. This led to a failure rate of **6%**. In comparison, 8 samples failed with Myriad MyChoice DX, representing a failure rate of 11%. Eliminating those failed samples, 63 left for the comparison.

Table 1: Distribution per center of samples and results used for the comparison study.

Center	Sequenced samples	Published samples (good QC)	Not Valuable from Myriad	Samples usable for the correlation analysis	Concordant HRD	Conflicting HRD
Hospital U. 12 de Octubre	18	18	3	15	14	1
Hospital Vall d'Hebron	19	17	2	15	15	0
Fundación Jiménez Díaz	22	19	2	17	14	3
Hospital Clínico U. Santiago	17	17	1	16	13	3
<b>TOTAL</b>	<b>76</b>	<b>71</b>	<b>8</b>	<b>63</b>	<b>56</b>	<b>7</b>




Figure 5: Distribution of negative and positive (BRCA1, BRCA2 or HRR altered or HRR wild-type) HRD samples.

- Of the 63 samples (Figure 5), 31 (48.4%) were **HRD positive** and 32 (51.6%) were **HRD negative**. Of the 31 HRD positives, 8 had alteration in BRCA1 (13%), 6 in BRCA2 (9%), 6 in other HRR genes (9%) and 17% without HRR variant, which is in line with literature. Overall, the most mutated gene was TP53, followed by HRR genes (BRCA1, BRCA2, CCNE1,...).

- 56 out of the 63 samples showed concordant results (88,8%) (Figure 6). The remaining 7 samples with conflicting results; 4 positives and 3 negatives and 5 fell in the grey zone (40+-5). All 14 BRCA1/2 alterations were detected (100%). Both expected concordance scores were passed and the HRD analysis showed a sensitivity of 90% and a specificity of 88%.




Figure 6: Distribution of samples showing concordant or conflicting HRD results compared to Myriad MyChoice DX (A) and distribution of the conflicting GS (B).




Figure 7: Additional information given by OncoDEEP in HRD positive patients highlighting alterations linked to shorter responses to PARP inhibitors.

- OncoDEEP panel is **more than just a HRD test (Figure 7)**. 3 HRD positive samples showed alterations possibly impacting PARP inhibitors response and allowing the identification of patients who **need closer monitoring**.

- Therapeutical options in 27/32 HRD negative patients (85%)(Figure 8): 2 patients eligible for immunotherapy and 25 patients harboring alterations matching to recruiting clinical trials




Figure 8: Therapeutical alternatives given by OncoDEEP in HRD negative patients showing patients eligible for immunotherapy (A) or directly enrollable to clinical trials through detected alterations (B).