



OncoDNA[®]
THE CANCER THERANOSTIC COMPANY

European User Group Meeting OncoDEEP[®] Kit


21st - 22nd October 2024



- Flexible and Innovative –
Delivering on the bigger picture

Koenraad Eycken

Being at the top of a process requires mastering every single step



Wetlab



Automation



Sequencing



Analysis



Clinical interpretation

We've designed a robust, adaptive workflow to fit most lab setups



DNA extraction impact

3938

Comprehensive Genomic Profiling of solid tumor patients with the OncoDEEP assay for broad analysis in clinical diagnostics

Hendrikus Jan Dubbink¹, Pieter-Jan Volders², Joni Van der Meulen³, Aaron De Cock³, Stefanie Vermeire⁴, Jacques Van Huysse⁵, Marie De Barys⁵, Gabriela Beniuga⁵, Zeliha Ozgur¹, Wendy de Leng⁶, Anne Jansen⁶, Ernst-Jan Speel⁷, Imke Demers⁷, Wilfred van Icken⁸, Brigitte Maes² and Guy Froyen^{2*}

¹ Dept. of Pathology, Erasmus MC Cancer Institute, Rotterdam, The Netherlands; ² Lab for Molecular Diagnostics, Jessa Hospital and University of Hasselt, Belgium; ³ Molecular Diagnostics Ghent University Hospital, Belgium; ⁴ Dept. of Pathology, AZ Sint-Jan Brugge AV, Bruges, Belgium; ⁵ Institute of Pathology and Genetics (IG), Gosselies, Belgium; ⁶ Dept. of Pathology, University Medical Centre Utrecht, The Netherlands; ⁷ Dept. of Pathology, Maastricht University Medical Center, The Netherlands; ⁸ Genomics Core Facility, Erasmus University Medical Center, Rotterdam, The Netherlands; * corresponding author.

Introduction

With the fast-growing number of recommended and required genomic biomarkers small gene panels have become vastly insufficient for most tumor types. Comprehensive Genomic Profiling (CGP) is amenable to screen for subtle nucleotide variants (SNVs and indels) in several hundred of cancer-related genes. Moreover, CGP can provide information on copy number variations (CNVs), gene fusions and tumor-agnostic genomic biomarkers including microsatellite instability (MSI), tumor mutation burden (TMB) and homologous recombination deficiency (HRD) for optimal clinical patient management with diagnostic, prognostic and therapeutic value in a wide variety of solid tumors. Only few CGP panels have been diagnostically validated in the clinic. Here, we report on an extensive multicentric (7 labs) comparative analysis of the novel CGP assay OncoDEEP from OncoDNA, with the diagnostically validated TSO500 assay from Illumina.

Table 1. Comparison of the number of genes for variant calling and the ability of biomarker detection for the TSO500 and OncoDEEP assays

Detection at DNA level	TSO500	OncoDEEP	in common
	# genes		
SNVs and indels	523	638	516
Amplifications	59	638	59
Deletions	0	xx	0
pan-tumor markers			
MSI	Yes	Yes	Yes
TMB	Yes	Yes	Yes
LOH	No	Yes	No
HRD	No*	Yes	No*
*separate HRD panel has to be added to the assay			
Detection at RNA level	TSO500	OncoDEEP	in common
	# genes		
Fusions	55	11	11
Splice variants	3	9	3

Materials and Methods

Both assays were performed as described in the user guides of both assays. In total, 234 diagnostic DNA and RNA samples with known TSO500 data were analyzed with the OncoDEEP assay. In addition, 12 and 8 reference DNA and RNA samples resp., were analysed for exon skipping and gene fusion detection by most laboratories. The diagnostic samples included many different tumor types, representative of the real life situation in the NGS diagnostic centers. Pooled libraries of both assays were sequenced on a NextSeq500/550 or NovaSeq instrument (Illumina). The major differences between both assays are listed in Table 2.

Table 2. Comparison of TSO500 and OncoDEEP assay features

	TSO500 (Illumina)	OncoDEEP (OncoDNA)
	Recommended input	DNA: 80 ng RNA: 40 ng
Pre-analytics		
DNA Fragmentation method	Shearing	Enzymatic
Normalisation	With beads	Quantification and dilution
Library prep		
Pooling before hybrid	No	Yes (8 samples)
# Hybridization steps	2	1
Hybridization capture		
Sequencing on a NextSeq550		
Read length	101 bp	74 bp
#Samples per run	8 (DNA + RNA)	24 (DNA + RNA)
Flowcell/NextSeq550Dx	High Output Kit v2.5 (300 Cycles)	High Output Kit v2.5 (150 Cycles)
Data analysis		
Secondary analysis	TSO500 local app	OncoKDM (OncoDNA)
Tertiary analysis	not available	OncoKDM (OncoDNA)
Hands-on-time and cost		
Hands-on-time	5 h	4 h
Cost/sample	\$\$\$	\$

Conclusions

- The OncoDEEP assay can efficiently detect somatic variants and CNV's in a broad range of solid tumor samples.
- Pan-cancer biomarker analysis (MSI, TMB, HRD) is mostly in line with TSO500 but threshold settings recommended.
- Gene fusion detection at the RNA level is only possible for 9 most common driver genes in solid tumors.
- The targeting capture provides a uniform selection of the targeted regions in a single hybridisation step.
- The assay includes variant classification and interpretation (OncoKDM) generating patient reports with local clinical trial selection.
- Successful validation (precision, sensitivity, specificity, accuracy, limit-of-detection) of the OncoDEEP assay was performed on 82 diagnostic and 8 reference samples for its implementation in clinical diagnostics.

Results

General comments OncoDEEP:

- The mean coverage of the samples is more uniform than for TSO500 (Figure 1) thereby allowing to pool 2- to 3-times more samples.
- The QC metrics are more prone to the DNA extraction procedure.

Variant detection with OncoDEEP:

- High concordance between both assays for SNP and indel detection (Figure 2). Missed variants (15%) were due to:
 - Insufficient coverage due to low QC sample metric (8%)
 - VAF did not reach the threshold of 5% (3%)
 - Too stringent filter settings e.g. >20 variant reads required (2%)
 - Classified as VUS instead of (Likely) Pathogenic (1%)
 - Different usage of NM_ transcript number (1%)

• Amplifications, if >6-fold change, were all concordant; LOH can not be assessed since it is not validated for TSO500.

• Gene fusions (82) or exon skipping (16) events were concordant in 78 cases (80%). In the undetected cases it was due to:

- Very low number of supporting reads in TSO500 (4 cases)
- Absence of driver gene in OncoDEEP panel (14 cases)
- Unknown reason (2 cases)

Pan-cancer biomarkers:

- In general good agreement between both assays
- Compared to TSO500, the values for MSI calling tend to be somewhat higher (Figure 3) while somewhat lower for TMB
- HRD could only be compared for few samples but were concordant




Figure 1. Distribution of the median coverage obtained with TSO500 (blue) and OncoDEEP (red) (SD 217 vs 145) showing the higher uniformity of the OncoDEEP capture.




Figure 2. Of the 678 TSO500 (Likely) Pathogenic SNV and indel variants 85% were also detected with the OncoDEEP assay (green). However, 102 (red) were absent due to a variety of reasons (see text).




Figure 3. MSI positive ratio plot (log-scale) for FPPE samples analyzed with TSO500 and OncoDEEP. Using a threshold of 20%, only 2 discordant calls (red) were present.




Figure 4. Comparison of TMB values obtained with TSO500 and OncoDEEP (log-scale). With a threshold of 16 mut/Mb 11 discordant calls (red) were found.

- Multicenter OncoDEEP – TSO500 comparison
- Data sets needed to be discarded due to bad extraction techniques

DNA extraction


DNA extraction can have a big impact on the final result.



AS1450		Aligned Reads	Reads Length (bp)	Duplicate (%)	On Target reads (%)	Uniformity (%)	GC Content (%)	Mean Coverage	Target Not Covered (%)	Target Low covered (%)	Mode Insert Size (bp)	Mean Insert Size (bp)
	KDM43635	24584834	99	25	40,2	83	51,6	297	0,02	4,40	138	198
	KDM43637	17095704	97	27	46,5	97	48,8	254	0,02	1,90	114	159
	KDM43638	18847627	94	21	47,0	98	50,5	279	0,02	1,20	101	131
	KDM43639	31309418	96	22	46,4	70	52,8	450	0,01	5,00	115	150
	KDM43640	19740875	97	21	48,7	98	48,9	308	0,02	0,80	113	156
	KDM43641	34433258	96	21	45,2	67	53,2	475	0,02	3,20	113	154
	KDM43642	25718423	98	21	45,4	97	49,8	365	0,02	0,40	142	183

AS1720	RD-23-112-Or	Aligned Reads	Reads Length (bp)	Duplicate (%)	On Target reads (%)	Uniformity (%)	GC Content (%)	Mean Coverage	Target Not Covered (%)	Target Low covered (%)	Mode Insert Size (bp)	Mean Insert Size (bp)
	KDM43889	26370315	98	18	37,6	99	49,1	320	0,02	0,50	115	171
	KDM43888	22360048	97	16	39,6	99	49,3	291	0,03	0,70	114	143
	KDM43887	26911421	97	18	38,3	99	49,2	338	0,02	0,50	114	151
	KDM43886	31139055	98	17	38,4	99	49,5	388	0,03	0,30	120	163
	KDM43885	25986585	98	17	38,7	99	48,8	328	0,02	0,50	122	162
	KDM43884	21314142	97	17	38,1	99	48,9	266	0,03	1,10	110	146
	KDM43883	25971479	98	18	38,0	99	49,1	320	0,02	0,50	122	164

Automation




Customized support

- Protocol review with the company
- Proposals of workflow fitting our assay requirements but getting around the tricky steps to automate
- Active role in the script development and testing phases



Sequencers



illumina®

MGI

Element
Biosciences

What has been done

- Clinical samples, carrying known variants/fusion/genomic signatures previously detected with OncoDEEP® kit combined with Illumina sequencing, were selected (**Table 1**). Reference standards were also used: 6 from Horizon for variants & fusions and 3 from Seracare for tumor mutational burden (TMB). DNA inputs of 6ng, 50ng & 100ng were used and some samples were replicated. In total, 90 samples (54 clinical & 36 reference samples) were sequenced to evaluate the performance of MGI System.

Table 1: Samples used to compare OncoDEEP® panel based on Illumina or MGI sequencer

Samples	Aim	n
Clinical samples	Variant/fusion calling and MSI	34
Ovarian clinical samples*	HRD	20
HD678, HD789, HD798, HD799, HD803, HD832, GM24149	Variant/fusion calling	29
Seracare samples	TMB	7

* FFPE ovarian clinical samples coming from the Hospital 12 Octubre (Madrid, Spain)

Gene	Variant (VAF%) ^a	HD798 (mild)				HD799 (moderate)				HD803 (severe)				HD832					
		Illumina		MGI		Illumina		MGI		Illumina		MGI		Illumina			MGI		
		50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	100ng	50ng **	6ng	100ng	50ng **	6ng
BRAF	V600E (10.5/10.7)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
KIT	D816V (10)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
EGFR	E746_A7 50 del (2/1.9)	✓	✓	✓	✓	✓	x [▲]	✓	✓	✓	x [▲]	✓	x [▲]	✓	✓	✓	✓	✓ [⊙]	✓ [⊙]
	L858R (3/2.8)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	T790M (1/0.9)	✓	✓ [■]	✓ [⊙]	✓ [■]	✓	x [▲]	✓	x [▲]	✓	x [▲]	✓	✓	✓ [■]	✓ [⊙]	✓	✓	✓ [⊙]	✓ [■]
	G719S (24.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
KRAS	G13D (15)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	G12D (6/6.3)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NRAS	Q61K (12.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
PIK3CA	H1047R (17.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	HD803	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

^a VAF in green is specific for Reference HD832

* Duplicated

** Quadruplicate

▲ Detected but discarded due to failing QC cut-off

■ Detected but not reported due to VAF <1%

⊙ 1 duplicate detected but not reported due to VAF <1%

⊙⊙ 2 duplicates detected but not reported due to VAF <1%

Element Bioscience

Gene	Variant (VAF%) [⊗]	HD798 (mild)				HD799 (moderate)				HD803 (severe)				HD832					
		Illumina		EB		Illumina		EB		Illumina		EB		Illumina			EB		
		50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	50ng *	6ng	100ng	50ng **	6ng	100ng	50ng **	6ng
BRAF	V600E (10.5/10.7)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
KIT	D816V (10)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
EGFR	E746_A7 50 del (2/1.9)	✓	✓	✓	✕ [▲]	✓	✕ [▲]	✓ [△]	✕ [▲]	✓	✕ [▲]	✓	✓	✓	✓	✓	✓ [■]	✓ ^{△△}	✕ [▲]
	L858R (3/2.8)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	T790M (1/0.9)	✓	✓ [■]	✓ [△]	✕ [▲]	✓	✕ [▲]	✓	✕ [▲]	✓	✕ [▲]	✓ [○]	✓	✓ [■]	✓ ^{○○}	✓	✓	✓	✓
	G719S (24.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
KRAS	G13D (15)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	G12D (6/6.3)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NRAS	Q61K (12.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
PIK3CA	H1047R (17.5)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	E545K (9/8.8)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

⊗ VAF in green is specific for Reference HD832

■ Detected but not reported due to VAF <1%

* Duplicated

○ 1 duplicate detected but not reported due to VAF <1%

** Quadruplicate

○○ 2 duplicates detected but not reported due to VAF <1%

▲ Detected but discarded due to failing QC cut-off

△ 1 duplicate detected but discarded due to failing QC cut-off

△△ 2 duplicates detected but discarded due to failing QC cut-off


Provider	Sequencer	FC type	Lanes	Total Reads/FC	Read Lengths	Data Output	Run Time	Nb OncoDEEP DNA+RNA samples ‡
Illumina	Next Seq 500 / 550	Mid Output	4	260M	PE75	16-19 Gb	~15 hrs	8
				260M	PE150	32-39 Gb	~26 hrs	8
		High Output	4	800M	PE75	50-60 Gb	~18 hrs	16
				800M	PE150	100-120 Gb	~29 hrs	24
	NextSeq 2000	P1	1	200 M	PE150	30 Gb	~19 hrs	8
		P2	1	800M	PE100	80 Gb	~21 hrs	24
					PE150	120 Gb	~29 hrs	24-32
	P3	1	2400M	PE100	240 Gb	~33 hrs	72-96	
				PE150	360 Gb	~48 hrs	96	
	Novaseq 6000	SP	2	1300-1600M	PE100	134-167 Gb	~19 hrs	48
		SP	2		PE150	200-250 Gb	~25 hrs	48-56
		S1	2	2600-3200M	PE100	266-333 Gb	~19 hrs	112
		S1	2		PE150	400-500 Gb	~25 hrs	112-120
		S2	2	6600-8200M	PE100	667-833 Gb	~25 hrs	232*
S2		2	PE150		1000-1250 Gb	~36 hrs	240*	
Novaseq X	1.5B	2	3200M	PE100	330 Gb	~18 hrs	112	
	1.5B	2		PE150	500 Gb	~21 hrs	112-120	
MGI	DNBSEQ-G400	FCS	2	1100M	PE100	110 Gb	~22/26 hrs	16-24
					PE150	165 Gb	~31/37 hrs	24*
		FCL	4	3000M~3600M	PE100	300-360 Gb	~35/38 hrs	96
					PE150	450-540 Gb	~50/56 hrs	96-104*
	DNBSEQ-T7	DNBSEQ-T7RS HT V3.0	1	11600M	PE100	1160 Gb~4640 Gb	16hrs ~18hrs	368*
					PE150	1750 Gb~7000 Gb	22hrs ~24hrs	368-376*
		DNBSEQ-T7RS HT V2.0	1	10000M	PE100	1000 Gb~4000 Gb	20hrs ~22hrs	320*
					PE150	1500 Gb~6000 Gb	24hrs ~30hrs	320-328*
DNBSEQ-G50	FCL	1	1000M	PE100	~100 Gb	~26 hrs	24*	
				PE150	~150 Gb	~40 hrs	24-32*	
DNBSEQ-G99	FC	1	160M	PE150	24-48 Gb	12 hrs	4-5*	
Element	AVITI	Low Output	2	500M	PE150	75 Gb	27 hrs	8*
		Medium Output	2	1000M	PE75	75 Gb	20 hrs	24
					PE150	150 Gb	31 hrs	32*
		High Output	2	2000M	PE75	150 Gb	24 hrs	48
PE150	300 Gb				38 hrs	56-64*		
PacBio	Onso			800-1000M	PE100	80-100 Gb	32 hrs	24**
					PE150	120-150 Gb	48 hrs	24-32**

‡ Sample number includes DNA samples only or DNA samples plus the same number of RNA samples (eg; 8 would mean 8 DNA samples or 8 DNA samples plus 8 RNA samples)

*Please contact the support team to discuss on the feasible number of samples

** Under evaluation

Versatile



Solid Biopsy



Solid Tumor WES

RNA Panels

LargeRNA

FocusRNA

Liquid Biopsy



Monitoring



+ New developments

Most comprehensive genomic profiling panel

Designed by oncology experts, the OncoDEEP® panel contains **the most relevant and complete cancer gene panel**. Over time this panel was optimized to include all clinically-relevant oncology targets. The DNA and RNA panels comprise **638 and 22 genes** reporting **genomic alterations** (SNV, insertion, deletion, CNV, fusions, unusual splicing, ...) and **complex genomic signature** (HRD, MSI and TMB).



Article

[HRD assessment for clinical management using the OncoDEEP® Kit](#)

DNA (638 GENES)

Genes alterations

SNV (including TERT PROMOTER)

CNV

DELETION

INSERTION

LOH

Complex signatures

HRD

MSI

TMB

Translocation (intronic)

ALK

ROS1

RET

MET-EX14

RNA (22 DRIVER GENES)

Gene fusions

ALK

ROS1

RET

FGFR1

FGFR2

FGFR3

NTRK1

NTRK2

NTRK3

BRAF

NRG1

EWSR1

TPR52

Unusual Splicing

BRCA1

BRCA2

PTEN

AR

EGFR

ERBB2

MET


PALB2

RB1

Targeted Liquid Biopsy to Support Your Clinical Decisions

Designed by oncology experts, the **OncoSELECT** panel contains the most relevant cancer gene panel associated with metastatic solid cancer, including **all clinically-relevant oncology targets**.

The panel is composed of **74 genes**, reporting genomic alterations (SNV, insertion, deletion, intron aberrations such as gene translocations and unusual splicing) and TERT promoter



Targeted Liquid Biopsy to Support Your Clinical Decisions

Optimized process to achieve **greater sensitivity and specificity**

VOF =0,25%

VOF =0,1%

98,5 > % sensitivity

87,5% > % sensitivity

99,9 > % specificity

99,9 > % specificity

99,9 > % accuracy

99,9 > % accuracy

Translocation - 24 genes

*** All CDS - 28 genes**

AKT1	BRCA2*	EML4	FGFR3*	IDH2	NBN*	PIK3CA	SLC34A2
ALK	BRIP1*	ERBB2	FOXL2	KEAP1*	NF1*	POLE	STK11*
AR	CD74	ERBB4	GNA11	KIF5B	NPM1	PTEN*	pTERT
ARAF	CDK12*	ESR1	GNAQ	KIT	NRAS	RAD51B*	TP53
ARID1A*	CDKN2A*	ETV6	GNAS	KRAS	NRG1	RAD51C*	
ATM*	CHEK1*	EZR	H3F3A	MAP2K1	NTRK1	RAD51D*	
ATR*	CHEK2*	FANCA*	H3F3B	MET	NTRK2	RAD54L*	
BARD1*	CTNNB1	FANCL*	H3C2	MLH1*	NTRK3	RET	
BRAF	DICER1	FGFR1*	HRAS	MRE11A*	PALB2*	ROS1	
BRCA1*	EGFR	FGFR2*	IDH1	MYOD1	PDGFRA	SDC4	



**** Detection of MSI and CNV****

Details


- Gene list: *AKT1* (3), ***ALK***(20-25), ***BRAF*** (6-8, 11, 12, 14, 15), *BRCA1**, *BRCA2**, ***CD74***, *CDKN2A**, *CTNNB1* (3), *DICER1* (24,25), ***EGFR*** (18-21), ***EML4***, *ERBB2* (3, 6-8, 12, 16-21), *ERBB4* (10, 12), *ESR1* (5-8), ***ETV6***, ***EZR***, ***FGFR1****, ***FGFR2****, ***FGFR3****, *FOXL2* (1), *GNA11* (4,5), *GNAQ* (4,5), *GNAS* (8), *H3-3A* (*H3F3A*) (2), *H3-3B* (*H3F3B*) (2), *H3C2* (*HIST1H3B*) (1), *HRAS* (2-5), *IDH1* (4), *IDH2* (4), *KEAP1**, ***KIF5B***, *KIT* (8, 9, 11, 13, 14, 17-19), *KRAS* (2-4), *MAP2K1* (2,3), ***MET*** (13-20), *MYOD1* (1), ***NPM1***, *NRAS* (2-4), ***NRG1***, ***NTRK1*** (14, 15), ***NTRK2*** (12, 13), ***NTRK3*** (16, 17), *PDGFRA* (12, 14, 18), *PIK3CA* (2, 3, 5-8, 10, 14, 19, 21), *POLE**, *PTEN**, ***RET*** (8, 10-16), ***ROS1*** (37, 38, 40), ***SDC4***, ***SLC34A2***, *STK11**, *TERT* (promotor), *TP53**
 - *AR*, *ARAF*, *ARID1A**, *ATM**, *ATR**, *BARD1**, *BRIP1**, *CDK12**, *CHEK1**, *CHEK2**, *FANCA**, *FANCL**, *MLH1**, *MRE11A**, *NBN**, *NF1**, *PALB2**, *RAD51B**, *RAD51C**, *RAD51D**, *RAD54L**
-
- MSI, genomic backbone for CNV detection (amplification, LOH, deletion...)

*Whole gene


Translocations

(X) exon targeted for SNV/Indels


Use cases – Start combing



Start combining – scenario 1 – low throughput



Start combining – scenario 2 – Medium throughput




(800 samples / year)




+



(300 samples / year)



illumina®




NextSeq 550/550 high 2x75/2x150 for OncoDEEP
 NextSeq 500/550 high 2x150 for OncoSelect
 NextSeq 1000/2000 P2 2x100/2x150 for OncoDEEP
 NextSeq 1000/2000 P2 2x150 for OncoSELECT
 NextSeq 2000 P3 2x150 combining OncoDEEP and OncoSELECT

50 flow cells / year
 50 flow cells / year
 50 flow cells / year
 50 flow cells / year
 50 flow cells / year

OR

Element Biosciences




EB Aviti Medium output 2x75 for OncoDEEP
 EB Aviti Medium output 2x150 for OncoSELECT
 EB Aviti High output 2x150 combining OncoDEEP and OncoSELECT

50 flow cells / year
 50 flow cells / year
 50 flow cells / year

OR


MGI



DBSNQ400 MGI FCS 2x100 for OncoDEEP
 DBSNQ400 MGI FCS 2x150 for OncoSELECT
 DBSNQ400 MGI FCL 2x150 combining OncoDEEP and OncoSELECT

50 flow cells / year
 50 flow cells / year
 50 flow cells / year

Start combining – scenario 3 – high throughput



(1200 samples / year)




DNA


+



(600 samples / year)



illumina®



NextSeq 500/550 High 2x75/2x150 for OncoDEEP

2x NextSeq 550 high for OncoSelect

NextSeq 1000/2000 P2 2x100/2x150 for OncoDEEP

2x NextSeq 1000/2000 P2 2x150 for OncoSELECT

NextSeq 2000 P3 2x150 combining OncoDEEP and OncoSELECT

50 flow cells / year

100 flow cells / year


50 flow cells / year

100 flow cells / year

50 flow cells / year

OR

Element Biosciences



EB Aviti Medium output 2x75/2x150 for OncoDEEP

EB Aviti High output 2x150 for OncoSELECT


50 flow cells / year

50 flow cells / year

50 flow cells / year

OR

MGI



DBSNQ400 MGI FCS 2x150 for OncoDEEP


DBSNQ400 MGI FCL 2x150 combining OncoDEEP and OncoSELECT

50 flow cells / year


50 flow cells / year

50 flow cells / year

Our workflow offers multiple combinations to optimize efficiency and cost...



Even when the backlog become critical



Use Case – High Throughput lab

- **The ask**

- Medium size targeted tumor panel (SNV, InDel, CNV, fusion, MSI, PgX) and small RNA panel 2000 samples
- Large CGP tumor panel (SNV, InDel, CNV, fusion, MSI, PgX; HRD) and large RNA panel 800
- Targeted tumor panel for hematological samples (SNV, InDel, fusion) 672
- CF DNA panel 500 samples

- **The OncoDNA proposal**


- OncoDEEP for all solid tumor panels and hematological samples
- OncoSelect and potentially OncoDEEP on Liquid Biopsy panel
- Large RNA panel

		OncoDEEP Library Preparation and Capture - up to 48 samples batch	Duration	Hands-on Time Manual	Hands-on Time Automated
D1 8h30	STEP 1	DNA Fragmentation, End-Repair and dA-tailing dA-tailed DNA samples	1 hour	30 min.	40 min.
	STEP 2	Universal Adapter Ligation and Clean-up Ligation products	1 hour	15 min. 30 min. Clean-up	
	STEP 3	PCR Amplification with UDI Primers, Clean-up and QC Individual indexed libraries	1 hour	15 min. 30 min. Clean-up 15 min. QC	
	STEP 4	Preparation of libraries and pools for hybridization Pool of Individual indexed libraries	1 hour	15 min. + 15 min.	20 min.
	STEP 5	Capture probes hybridization Hybridized targets in solution	16 hour	15 min.	
D1 15h30					
D2 8h30	STEP 6	Hybridized targets binding to the streptavidin beads Captured targets on beads	1,5 hour	30 min. 30 min. Clean-up	15 min.
	STEP 7	Post-capture PCR amplification, clean-up and QC Enriched libraries ready to be sequenced	1 hour	30 min.	
D2 14H					

Why is this important

		Monday	Tuesday	Wednesday	Thursday	Friday
Automation	AM		Lib Prep RNA	Capture Washes & Post-Capture Amplification	Lib Prep RNA	Capture Washes & Post-Capture Amplification
	PM	Lib Prep DNA	Hybridization and Capture Launch	Lib Prep DNA	Hybridization and Capture Launch	
Manual		cDNA synthesis		Sequencing Run launch		Sequencing Run launch
				cDNA synthesis		

We are your partner to navigate a complex ecosystem



Data Data Data – Live demo

